


bsm-62728R**[Primary Antibody]****Centaurin beta 2 Recombinant Rabbit mAb****BioSS**
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —**Host:** Rabbit**Isotype:** IgG**Clonality:** Recombinant**GeneID:** 23527**SWISS:** Q15057**Target:** Centaurin beta 2**Immunogen:** A synthesized peptide derived from human ACAP2: 650-715.**Purification:** affinity purified by Protein A**Storage:** 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol.
Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.**Background:** GTPase-activating protein (GAP) for ADP ribosylation factor 6 (ARF6).**Applications:** WB (1:1000-2000)**Reactivity:** Human, Mouse, Rat**Predicted
MW.:** 88**Subcellular
Location:** Cytoplasm**— VALIDATION IMAGES —**

Western blot analysis of 293 cell lysate. Using Centaurin beta 2 (bsm-62728R) monoclonal antibody at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.