

bsm-62747R**[Primary Antibody]****Bioss**
ANTIBODIES

www.bioss.com.cn


sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

PWP1 Recombinant Rabbit mAb**— DATASHEET —**

Host: Rabbit	Isotype: IgG	Applications: WB (1:1000-2000)
Clonality: Recombinant		Reactivity: Human, Mouse, Rat
GeneID: 11137	SWISS: Q13610	
Target: PWP1		
Immunogen: A synthesized peptide derived from human PWP1: 460-501.		Predicted MW.: 56
Purification: affinity purified by Protein A		
Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.		
Background: Chromatin-associated factor that regulates transcription.Regulates Pol I-mediated rRNA biogenesis and, probably, Pol III-mediated transcription.Regulates the epigenetic status of rDNA.		

— VALIDATION IMAGES —

Western blot analysis of HepG2 cell lysate. Using PWP1 (bsm-62747R) monoclonal antibody at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.