


**bsm-62747R****[ Primary Antibody ]****Bioss**  
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

**PWP1 Recombinant Rabbit mAb****— DATASHEET —****Host:** Rabbit**Isotype:** IgG**Clonality:** Recombinant**GeneID:** 11137**SWISS:** Q13610**Target:** PWP1**Immunogen:** A synthesized peptide derived from human PWP1: 460-501.**Purification:** affinity purified by Protein A**Storage:** 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol.  
Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.**Background:** Chromatin-associated factor that regulates transcription.Regulates Pol I-mediated rRNA biogenesis and, probably, Pol III-mediated transcription.Regulates the epigenetic status of rDNA.**Applications:** WB (1:1000-1:2000)**Reactivity:** Human, Mouse, Rat**Predicted MW.:** 56**— VALIDATION IMAGES —**

Western blot analysis of HepG2 cell lysate. Using PWP1 (bsm-62747R) monoclonal antibody at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.