

bsm-62849R

[Primary Antibody]

FRG1 Recombinant Rabbit mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:1000-2000)
Clonality: Recombinant		Reactivity: Human, Mouse, Rat
GeneID: 2483	SWISS: Q14331	
Target: FRG1		
Immunogen: A synthesized peptide derived from human FRG1: 174-226.		Predicted MW.: 29
Purification: affinity purified by Protein A		Subcellular Location: Nucleus
Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.		
Background: May have a role in processing of pre-rRNA or in the assembly of rRNA into ribosomal subunits. May be involved in pre-mRNA splicing.		