

bsm-62849R

[Primary Antibody]

FRG1 Recombinant Rabbit mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

| | | |
|--|----------------------|---|
| Host: Rabbit | Isotype: IgG | Applications: WB (1:1000-1:2000) |
| Clonality: Recombinant | | Reactivity: Human, Mouse, Rat |
| GeneID: 2483 | SWISS: Q14331 | |
| Target: FRG1 | | |
| Immunogen: A synthesized peptide derived from human FRG1: 174-226. | | Predicted MW.: 29 |
| Purification: affinity purified by Protein A | | Subcellular Location: Nucleus |
| Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles. | | |
| Background: May have a role in processing of pre-rRNA or in the assembly of rRNA into ribosomal subunits. May be involved in pre-mRNA splicing. | | |