

**bsm-62854R**

**[ Primary Antibody ]**

## PTN Recombinant Rabbit mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

**Host:** Rabbit

**Isotype:** IgG

**Clonality:** Recombinant

**GeneID:** 5764

**SWISS:** P21246

**Target:** PTN

**Immunogen:** A synthesized peptide derived from human PTN: 110-168/168.

**Purification:** affinity purified by Protein A

**Storage:** 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol.  
Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.

**Background:** Heparin binding mitogenic protein. Has neurite extension activity.

**Applications:** WB (1:500-2000)

**Reactivity:** Human, Mouse, Rat

**Predicted  
MW.:** 19

**Subcellular  
Location:** Secreted