

bsm-62878R

[Primary Antibody]

NELFB Recombinant Rabbit mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Clonality: Recombinant

GeneID: 25920

SWISS: Q8WX92

Target: NELFB

Immunogen: A synthesized peptide derived from human NELF B: 1-200/580.

Purification: affinity purified by Protein A

Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.

Background: Essential component of the NELF complex, a complex that negatively regulates the elongation of transcription by RNA polymerase II. The NELF complex, which acts via an association with the DSIF complex and causes transcriptional pausing, is counteracted by the P-TEFb kinase complex.

Applications: WB (1:1000-1:2000)

IHC-P (1:100-1:200)

IHC-F (1:100-1:200)

IF (1:50-1:200)

ICC/IF (1:50-1:200)

Reactivity: Human

Predicted MW.: 66