

**bsm-62882R**

**[ Primary Antibody ]**

## SART1 Recombinant Rabbit mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

**Host:** Rabbit

**Isotype:** IgG

**Clonality:** Recombinant

**GeneID:** 9092

**SWISS:** O43290

**Target:** SART1

**Immunogen:** A synthesized peptide derived from human SART1: 660-735/800.

**Purification:** affinity purified by Protein A

**Storage:** 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.

**Background:** Plays a role in mRNA splicing as a component of the U4/U6-U5 tri-snRNP, one of the building blocks of the spliceosome. May also bind to DNA.

**Applications:** **WB** (1:1000-5000)

**Flow-Cyt** (1:50-100)

**ICC/IF** (1:50-200)

**Reactivity:** Human, Mouse, Rat

**Predicted  
MW.:** 90

**Subcellular  
Location:** Nucleus