

bsm-62882R**[Primary Antibody]**

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

SART1 Recombinant Rabbit mAb**— DATASHEET —****Host:** Rabbit**Isotype:** IgG**Clonality:** Recombinant**GeneID:** 9092**SWISS:** O43290**Target:** SART1**Immunogen:** A synthesized peptide derived from human SART1: 660-735/800.**Purification:** affinity purified by Protein A**Storage:** 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.**Background:** Plays a role in mRNA splicing as a component of the U4/U6-U5 tri-snRNP, one of the building blocks of the spliceosome. May also bind to DNA.**Applications:** **WB** (1:1000-1:5000)
Flow-Cyt (1:20-1:100)
ICC/IF (1:50-1:200)**Reactivity:** Human, Mouse, Rat**Predicted
MW.:** 90**Subcellular
Location:** Nucleus