

bsm-63024R

[Primary Antibody]

Thiophosphate ester Recombinant Rabbit mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Applications: **IP** (1:20-50)

Clonality: Recombinant

Target: Thiophosphate ester

Immunogen: KLH conjugated Thiophosphate ester.: Thiophosphate ester(硫代磷酸酯).

Purification: affinity purified by Protein A

Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol.
Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.

Background: Thiophosphate ester is used as fungicide.