

**bsm-63024R**

**[ Primary Antibody ]**

## Thiophosphate ester Recombinant Rabbit mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

**Host:** Rabbit

**Isotype:** IgG

**Clonality:** Recombinant

**Target:** Thiophosphate ester

**Immunogen:** A synthesized peptide derived from Thiophosphate ester:  
Thiophosphate ester(硫代磷酸酯).

**Purification:** affinity purified by Protein A

**Storage:** 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol.  
Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.

**Background:** Thiophosphate ester is used as fungicide.

**Applications:** IP (1:20-1:50)