
Human GFAP Recombinant Mouse mAb

Cat. No:	V5603
Name:	Human GFAP Recombinant Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Recombinant
Clone No:	3A2
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	(predicted: Human)
Format:	Liquid
Buffer:	PBS (pH7.4).
Applications:	ELISA
Antigen:	Recombinant Human GFAP Protein
Pairing Information:	V5607&V5608-V5603