

---

## Human GFAP Mouse mAb

Cat. No:	V5604
Name:	Human GFAP Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	7G5
Host/Source:	Mouse
Subtype:	IgG2b
Reactivity:	(predicted:&nbsp;Human)
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	PBS (pH7.4).
Applications:	ELISA
Antigen:	Recombinant Human GFAP Protein ( <a href="#">bs-42001P</a> )
Pairing Information:	
V5604-V5606	