

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn 400-901-9800

Human GFAP Recombinant Mouse mAb

Cat. No: V5604

Name: Human GFAP Recombinant Mouse mAb

Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Clone type: Recombinant

Clone No: 7G5
Host/Source: Mouse
Subtype: IgG2b

Reactivity: (predicted: Human)

Purity: >90% as determined by SDS-PAGE

Format: Liquid

Buffer: PBS (pH7.4).

Applications: ELISA

Antigen: Recombinant Human GFAP Protein (bs-42001P)

Pairing Information:

V5604-V5606