


bsm-33007R**[Primary Antibody]****Bioss**
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

GST tag Recombinant Rabbit mAb**— DATASHEET —****Host:** Rabbit**Clonality:** Recombinant**CloneNo.:** 1C14**Target:** GST tag**Purification:** affinity purified by Protein A**Storage:** Size : 50ul/100ul/500ul0.01M TBS(pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Shipped at 4°C. Store at -20 °C for one year. Avoid repeated
freeze/thaw cycles.**Background:** GST is a 26kDa protein encoded by the parasitic helminth
Schistosoma japonicum and widely used in the pGEX family of GST
plasmid expression vectors as a fusion protein with foreign proteins.**Applications:** **WB** (1:5000-100000)**ELISA** (1:1000-5000)**Reactivity:** Species independent**Subcellular
Location:** Cytoplasm**— VALIDATION IMAGES —**

100 ng rmS100A7-GST protein (bs-42023P) per
lane probed with GST tag monoclonal antibody
respectively, unconjugated (bsm-33007R) at
1:50000 dilution and 4°C overnight incubation.
Followed by corresponding conjugated
secondary antibody incubation at r.t. for 60 min.