

## GST tag Recombinant Rabbit mAb

Catalog Number: bsm-33007R

Target Protein: GST tag

Form: Size : 50ul/100ul/500ul

Liquid

Host: Rabbit

Clonality: Recombinant

Clone No.: 1C14

Applications: WB (1:5000-100000)

Reactivity: Mouse, Species independent

Purification: affinity purified by Protein A


Storage: Size : 50ul/100ul/500ul

0.01M TBS(pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

**Background:** GST is a 26kDa protein encoded by the parasitic helminth *Schistosoma japonicum* and widely used in the pGEX family of GST plasmid expression vectors as a fusion protein with foreign proteins.

### VALIDATION IMAGES



100 ng rmS100A7-GST protein (bs-42023P) per lane probed with GST tag monoclonal antibody respectively, unconjugated (bsm-33007R) at 1:50000 dilution and 4°C overnight incubation. Followed by corresponding conjugated secondary antibody incubation at r.t. for 60 min.