

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn 400-901-9800

Bovine PAG1 Mouse mAb

Cat. No: V8302

Name: Bovine PAG1 Mouse mAb

Storage: Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Clone type: Monoclonal

Clone No: 9D10
Host/Source: Mouse
Subtype: IgG2a
Reactivity: Bovine

Purity: >90% as determined by SDS-PAGE

Concentration: 0.01M PBS(pH7.4)

Applications: ELISA

Antigen: Recombinant Bovine PAG1 Protein (bs-43189P)

Pairing Information: V8302-V8301&V8303