
Bovine PAG1 Mouse mAb

Cat. No:	V8302
Name:	Bovine PAG1 Mouse mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	9D10
Host/Source:	Mouse
Reactivity:	Bovine
Purity:	>90% as determined by SDS-PAGE
Concentration:	0.01M PBS(pH7.4)
Applications:	ELISA
Antigen:	Recombinant Bovine PAG1 Protein (bs-43189P)
Pairing Information:	
V8302-V8301&V8303	