

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn 400-901-9800

Bovine PAG1 Mouse mAb

Cat. No: V8303

Name: Bovine PAG1 Mouse mAb

Storage: Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Clone type: Monoclonal

Clone No: 12G5
Host/Source: Mouse
Subtype: IgG2a

Reactivity: (predicted: Bovine)

Purity: >90% as determined by SDS-PAGE

Format: Liquid

Buffer: 0.01M PBS(pH7.4)

Applications: ELISA, ICA

Antigen: Recombinant Bovine PAG1 Protein

Pairing Information:

V8302-V8303