
Human Serpin B3 Mouse mAb

Cat. No:	V7701
Name:	Human Serpin B3 Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	12D7
Host/Source:	Mouse
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Applications:	ELISA
Antigen:	Recombinant Human Serpin B3 Protein (bs-43190P)
Pairing Information:	
V7701-V7702	