

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn 400-901-9800

## **Human Serpin B3 Mouse mAb**

Cat. No: V7701

Name: Human Serpin B3 Mouse mAb

Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Clone type: Monoclonal

Clone No: 12D7
Host/Source: Mouse
Reactivity: Human

Purity: >90% as determined by SDS-PAGE

Format: Liquid

Buffer: PBS(pH=7.4)

Applications: ELISA

Antigen: Recombinant Human Serpin B3 Protein (bs-43190P)

Pairing Information:

V7701-V7702