

---

## Human Serpin B3 Mouse mAb

Cat. No:	V7704
Name:	Human Serpin B3 Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	2H6
Host/Source:	Mouse
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	PBS (pH=7.4)
Applications:	ICA
Antigen:	Recombinant Human Serpin B3 Protein ( <a href="#">bs-43190P</a> )
Pairing Information:	
V7703-V7704	