
Human Serpin B3 Mouse mAb

| | |
|----------------------|---|
| Cat. No: | V7705 |
| Name: | Human Serpin B3 Mouse mAb |
| Storage: | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| Clone type: | Monoclonal |
| Clone No: | 2E9 |
| Host/Source: | Mouse |
| Reactivity: | Human |
| Purity: | >90% as determined by SDS-PAGE |
| Format: | Liquid |
| Buffer: | PBS (pH=7.4) |
| Applications: | ICA |
| Antigen: | Recombinant Human Serpin B3 Protein (bs-43190P) |
| Pairing Information: | |
| V7705-V7706 | |