
Human Serpin B4 Mouse mAb

Cat. No:	V7603
Name:	Human Serpin B4 Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	7B7
Host/Source:	Mouse
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	PBS (pH=7.4)
Applications:	ELISA
Antigen:	Recombinant human Serpin B4 protein (bs-43150P)
Pairing Information:	
V7604-V7603	