
Dog Cys-C Mouse mAb

Cat. No:	V8101
Name:	Dog Cys-C Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	1F6
Host/Source:	Mouse
Subtype:	IgG2a
Reactivity:	Dog
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	PBS (pH=7.4)
Applications:	ELISA, ICA
Antigen:	Recombinant Dog Cys-C Protein (bs-43166P)
Pairing Information:	
V8101-V8102&V8103	