
Human SHBG Mouse mAb

Cat. No:	V6301
Name:	Human SHBG Mouse mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	6B6
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	PBS (pH=7.4)
Applications:	ELISA
Antigen:	Recombinant human SHBG protein (bs-43202P)
Pairing Information:	
V6301-V6302& V6303	