
Human HCG-beta Mouse mAb

Cat. No:	V2708
Name:	Human HCG-beta Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	3F4
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Applications:	ELISA, ICA
Antigen:	Recombinant Human HCG-beta Protein (bs-43172P)
Pairing Information:	
V2708-V2709&V2710	