

---

## Human HCG-beta Mouse mAb

Cat. No:	V2709
Name:	Human HCG-beta Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	1E3
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	PBS (pH=7.4)
Applications:	ELISA
Antigen:	Recombinant Human HCG-beta Protein ( <a href="#">bs-43172P</a> )
Pairing Information:	
V2707-V2709	