
Human HCG-beta Mouse mAb

Cat. No:	V2710
Name:	Human HCG-beta Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	4D3
Host/Source:	Mouse
Subtype:	IgG2b
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	PBS (pH=7.4)
Applications:	ELISA, ICA
Antigen:	Recombinant Human HCG-beta Protein (bs-43127P)
Pairing Information:	
V2708-V2710	