
Cat Cys-C Mouse mAb

Cat. No:	V8204
Name:	Cat Cys-C Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	5C7
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	(predicted: Cat)
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	PBS (pH=7.4)
Applications:	ELISA, ICA
Antigen:	Recombinant Cat Cys-C Protein (bs-43167P)
Pairing Information:	
V8202-V8204	