
Recombinant human Cadherin-1 protein, N-His

Catalog Number: bs-42419P

Concentration: >1mg/ml

Species: Human

AA Seq: 155-367/882

Predicted MW: 27.1 kDa

Tags: N-His

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC

Form: Liquid

Storage: 20mM Tris-HCl (pH=8.0)

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: This gene encodes a classical cadherin of the cadherin superfamily. Alternative splicing results in multiple transcript variants, at least one of which encodes a preproprotein that is proteolytically processed to generate the mature glycoprotein. This calcium-dependent cell-cell adhesion protein is comprised of five extracellular cadherin repeats, a transmembrane region and a highly conserved cytoplasmic tail. Mutations in this gene are correlated with gastric, breast, colorectal, thyroid and ovarian cancer. Loss of function of this gene is thought to contribute to cancer progression by increasing proliferation, invasion, and/or metastasis. The ectodomain of this protein mediates bacterial adhesion to mammalian cells and the cytoplasmic domain is required for internalization. This gene is present in a gene cluster with other members of the cadherin family on chromosome 16. [provided by RefSeq, Nov 2015]

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.