
Cortisol Mouse mAb

Cat. No:	V13402
Name:	Cortisol Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	4C2
Host/Source:	Mouse
Subtype:	IgG1
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	0.01M PBS(pH7.4).
Applications:	ELISA, ICA
Antigen:	Cortisol-KLH (bs-25463P)
Pairing Information:	
bs-25463PA-V13402	