
AFP Mouse mAb, HRP conjugated

Cat. No:	V2004-HRP
Name:	AFP Mouse mAb, HRP conjugated
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	1B1
Host/Source:	Mouse
Subtype:	IgG2a
Reactivity:	(predicted: Human)
Format:	Liquid
Concentration:	1mg/ml
Buffer:	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Applications:	ELISA
Antigen:	Recombinant Human AFP Protein