

---

## phospho-Tau (Thr217) Mouse mAb

Cat. No:	V5228
Name:	phospho-Tau (Thr217) Mouse mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	1G9
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	0.01M PBS(pH7.4)
Applications:	ELISA
Antigen:	KLH conjugated Synthesised phosphopeptide derived from human Tau around the phosphorylation site of Thr217 ( <a href="#">bs-2842P</a> )
Pairing Information:	
V5234&V5212-V5228	