
phospho-Tau (Thr217) Mouse mAb

Cat. No:	V5229
Name:	phospho-Tau (Thr217) Mouse mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	2B7
Host/Source:	Mouse
Subtype:	IgG2b
Purity:	Purified by Protein A
Format:	Liquid
Buffer:	0.01M PBS(pH7.4)
Antigen:	KLH conjugated Synthesised phosphopeptide derived from human Tau around the phosphorylation site of Thr217 (bs-2842P)
Pairing Information:	
V5234&V5212-V5229	