
Human GFAP Mouse mAb

Cat. No:	V5605
Name:	Human GFAP Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	8H10
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	PBS (pH7.4)
Applications:	ELISA
Antigen:	Recombinant Human GFAP Protein (bs-42001P)
Pairing Information:	
V5605-V5606	