

---

## Human GFAP Mouse mAb

|                      |                                                                                 |
|----------------------|---------------------------------------------------------------------------------|
| Cat. No:             | V5606                                                                           |
| Name:                | Human GFAP Mouse mAb                                                            |
| Storage:             | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| Clone type:          | Monoclonal                                                                      |
| Clone No:            | 12C12                                                                           |
| Host/Source:         | Mouse                                                                           |
| Subtype:             | IgG1                                                                            |
| Reactivity:          | Human                                                                           |
| Purity:              | >90% as determined by SDS-PAGE                                                  |
| Applications:        | ELISA                                                                           |
| Antigen:             | Recombinant Human GFAP Protein ( <a href="#">bs-42001P</a> )                    |
| Pairing Information: |                                                                                 |
| V5604&V5605-V5606    |                                                                                 |