

**bsm-25842M**

**[ Primary Antibody ]**

## zearalenone Mouse mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

**Host:** Mouse

**Isotype:** IgG

**Applications:** ELISA (1:5000-10000)

**Clonality:** Monoclonal

**Reactivity:** (predicted: Zearalenone)

**Target:** zearalenone

**Purification:** affinity purified by Protein A

**Storage:** 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.  
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

**Background:** Zearalenone is a mycotoxin affecting cereal grains for human and animal consumption.