
phospho-Tau (T212) Mouse mAb

Cat. No:	V5242
Name:	phospho-Tau (T212) Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	6E9
Host/Source:	Mouse
Subtype:	IgG2a
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	0.01M PBS(pH7.4)
Applications:	ELISA
Antigen:	KLH conjugated Synthesised phosphopeptide derived from human Tau around the phosphorylation site of Thr212 (bs-5419P)