

---

## Human UCHL1/PGP9.5 Recombinant Mouse mAb

Cat. No:	V5808
Name:	Human UCHL1/PGP9.5 Recombinant Mouse mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Recombinant
Clone No:	4C5
Host/Source:	Mouse
Subtype:	IgG2b
Reactivity:	(predicted:&nbsp;Human)
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	0.01M PBS(pH7.4)
Applications:	ELISA
Antigen:	Recombinant Human PGP9.5 Protein ( <a href="#">bs-41257P</a> )
Pairing Information:	V5807&V5809&V5810&V5811-V5808