
Human UCHL1/PGP9.5 Recombinant Rabbit mAb

Cat. No:	V5810
Name:	Human UCHL1/PGP9.5 Recombinant Rabbit mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Recombinant
Clone No:	13A5
Host/Source:	Rabbit
Subtype:	IgG
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	0.01M PBS(pH7.4)
Applications:	ELISA
Antigen:	Recombinant Human PGP9.5 Protein (bs-41257P)
Pairing Information:	
V5810-V5808	