
Vitamin B12 Mouse mAb

Cat. No:	V13701
Name:	Vitamin B12 Mouse mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	7B12
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Applications:	ELISA, ICA
Antigen:	VB12-KLH (D10788)
Pairing Information:	
bs-25860PB-V13701	