

---

## NF-L Mouse mAb

Cat. No:	V5512
Name:	NF-L Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	3E2
Host/Source:	Mouse
Subtype:	IgG2a
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA
Antigen:	Recombinant Human NF-L Protein ( <a href="#">bs-41217P</a> )
Pairing Information:	
V5508-V5512	