

bs-0464R

[Primary Antibody]

Hydroxyethyk starch Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Applications: ELISA (1:5000-10000)

Clonality: Polyclonal

Reactivity: (predicted: Hydroxyethyk starch)

Target: Hydroxyethyk starch

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Predicted MW.: 130 kDa

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.