

---

## Human CEA Mouse mAb, Biotin conjugated

Cat. No:	V1904-Bio
Name:	Human CEA Mouse mAb, Biotin conjugated
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	6C9
Host/Source:	Mouse
Subtype:	IgG
Reactivity:	(predicted:&nbsp;Human)
Format:	Liquid
Concentration:	1mg/ml
Buffer:	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Applications:	ELISA
Antigen:	KLH conjugated synthetic peptide derived from human Human CEA / Biotin ( <a href="#">bs-43141P</a> )