

---

## Rabbit Anti-Goat IgG H&L, BF594 conjugated

Catalog Number: bs-0294R-BF594

Target Protein: Rabbit Anti-Goat IgG H&L

Concentration: 2.0 mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IF (1:100-1000)

Reactivity: Goat

Purification: affinity purified by Protein A

Storage: 10 mM TBS (pH=7.4) with 1% BSA, 0.03% Proclin300 and 50% glycerol.

Store at -20°C for one year. Avoid repeated freeze/thaw cycles.