

Human VEGF165 Mouse mAb

Cat. No: V2603
 Name: Human VEGF165 Mouse mAb
 Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
 Clone type: Monoclonal
 Clone No: 3H7
 Host/Source: Mouse
 Subtype: IgG1
 Reactivity: (predicted: Human)
 Format: Liquid
 Buffer: 0.01M PBS (pH7.4).
 Applications: ELISA
 Antigen: Recombinant human VEGF165 Protein ([bs-41613P](#))
 Pairing Information:

ELISA		Detection	
		V2603	V2604
Capture	V2603	-	+++
	V2604	-	-