

bs-0239R

[Primary Antibody]

HIV gp41 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: ELISA (1:5000-10000)
Clonality: Polyclonal		Reactivity: (predicted: HIV)
Target: HIV gp41		
Immunogen: KLH conjugated synthetic peptide derived from Transmembrane protein gp41: 1-80/120.		Predicted MW.: 39/93 kDa
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		

— SELECTED CITATIONS —

- **[IF=4]** Wen-Qiang He. et al. IFI27 inhibits HIV-1 replication by degrading Gag protein through the ubiquitin-proteasome pathway. J VIROL. 2024 Oct 30 WB ;Human. 39475279