

bs-19033R

[Primary Antibody]

NBPF11 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) Reactivity: (predicted: Human) Predicted MW.: 99 kDa Subcellular Location: Cytoplasm
Clonality: Polyclonal		
GeneID: 200030	SWISS: Q86T75	
Target: NBPF11		
Immunogen: KLH conjugated synthetic peptide derived from human NBPF11: 781-865/865.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		