

**bs-4229R**

**[ Primary Antibody ]**

## **ZDHHC21 Rabbit pAb**



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

<b>Host:</b> Rabbit	<b>Isotype:</b> IgG	<b>Applications:</b> <b>WB</b> (1:500-2000) <b>ELISA</b> (1:5000-10000)
<b>Clonality:</b> Polyclonal		
<b>GeneID:</b> 340481	<b>SWISS:</b> Q8IVQ6	<b>Reactivity:</b> (predicted: Human, Mouse, Rat, Rabbit, Cow, Horse)
<b>Target:</b> ZDHHC21		
<b>Immunogen:</b> KLH conjugated synthetic peptide derived from human ZDHHC21: 131-240/265.		
<b>Purification:</b> affinity purified by Protein A		<b>Predicted MW.:</b> 31 kDa
<b>Concentration:</b> 1mg/ml		<b>Subcellular Location:</b> Cell membrane
<b>Storage:</b> 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
<b>Background:</b> ZDHHC21 belongs to the DHHC palmitoyltransferase family. The function of this protein has not been specifically defined.		